

# Certificate of Calibration TOC

**Model: PAT700**

**Serial Number: 1611**

**Calibration Date: 3/16/16**

**Calibration Due Date:** - *The PAT700 is designed to have long term stability. However, where required by the user, Anatel recommends re-calibration - two years after original calibration date or one year after installation date whichever is earlier.*

Anatel certifies that the above **Referenced Instrument** has been calibrated by comparison to or measurement against national standards as defined by the National Institute of Standards and Technology within the limitations of the Institute's Calibration Services, or by accepted values of natural physical constants, or by consensus standards or by ratio type measurements using self-calibrating techniques. The instrument is designed to comply with the latest TOC requirements specified by global Pharmacopeia including, but not limited to; USP United States Pharmacopeia, JP Japan Pharmacopeia, EP European Pharmacopeia, IP India Pharmacopeia, CP China Pharmacopeia, and KP Korean Pharmacopeia.

This instrument has been manufactured, inspected and tested in accordance with **ISO 9001:2008** Quality Standards. All records of inspection, testing and other supporting documentation are **Controlled Documents** and cannot be copied. However, these documents are available for on-site review. Hach's written description of its **Calibration Control System** is available upon request.

*Svet Thoms*  
\_\_\_\_\_  
Hach Technician Signature

*3/16/16*  
\_\_\_\_\_  
Date

**Anatel**<sup>®</sup>



**Be Right**<sup>™</sup>

Document #: 19002-00-PAT700-COC